

TC 900

Upper enclosure

For 900 mm-wide appliances without a steam condenser or drying unit.



- Includes lockable service door
- Wall panelling up to ceiling must be provided on site

EAN: 4002515991756 / Material number: 10880700 / Old Material Number: 69252003D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
Product category	
Cover	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in in. (mm)	33 1/8 (840)
External dimensions, net width in in. (mm)	33 5/16 (845)
External dimensions, net depth in in. (mm)	41 3/16 (1045)
External dimensions, gross height in in. (mm)	39 (990)
External dimensions, gross width in in. (mm)	33 5/16 (845)
External dimensions, gross depth in in. (mm)	47 1/4 (1200)
Net weight in lb (kg)	108 (48.90)
Gross weight in lb (kg)	141 (63.90)