

ProCare Shine 40 - 1 l

Rinse aid, 1 l

For rinsing dishes, cutlery and glasses.

EAN: 4002516054436 / Material number: 10923710 / Old Material Number: 61915101USA



- Slightly acidic, contains citric acid
- For optimum rinsing results
- For use in commercial dishwashers
- Liquid concentrate, recommendation: 0.02-0.05 % (0.2-0.5 ml/l)
- Specially developed for Miele Professional

Product category	
Fresh-water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Dispensing recommendation	
Dispensing recommendation per rinse	0.2-0.5 ml/l (0.02-0.05%)
pH value at 20°C (68°F)	1,9 - 3,0
Recommended water range in mmol/l	0-0.53
Product characteristics	
Medium type	Rinse aid
Aggregate form	Liquid concentrate
Ingredients	Sequestering agent (no
Suitable for glasses	•
Suitable for cutlery	•
Suitable for plastic	•
Suitable for stainless steel	•
Suitable for porcelain / ceramics	•
Suitable for aluminum	•
Suitable for silver	•
Gentle on decor	•
Prevents limescale and protects load items	•
Improves and shortens the drying process	•
Waste water neutral	•
Dimensions and weight	
External dimensions, net height in in. (mm)	10 9/16 (267)
External dimensions, net width in in. (mm)	3 1/16 (77)
External dimensions, net depth in in. (mm)	3 1/16 (77)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	3 3/16 (80)
External dimensions, gross depth in in. (mm)	3 3/16 (80)
Gross weight in lb (kg)	0.04 (1)
Fill volume in fl.oz. (ml)	33.81 (1000) fl oz (US) (Milliliter)
Compliance, test certificates, and symbols	
GHS / CLP compliant	•