

A 500

Mobile unit

With 1 level for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays
- Extendable to 2 levels with A 603

EAN: 4002516050421 / Material number: 10973880 / Old Material Number: 69650001D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Mobile unit for inserts	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Capacity	
DIN mesh trays [number]	3
SPRI mesh trays [number]	2
Dimensions and weight	
External dimensions, net height in in. (mm)	6 1/16 (153)
External dimensions, net width in in. (mm)	25 5/8 (650)
External dimensions, net depth in in. (mm)	31 1/8 (790)
External dimensions, gross height in in. (mm)	6 3/4 (170)
External dimensions, gross width in in. (mm)	26 (660)
External dimensions, gross depth in in. (mm)	31 15/16 (810)
Net weight in lb (kg)	11 (4.70)
Gross weight in lb (kg)	12 (5.00)