

A 501

Mobile unit

With 4 levels for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays
- Extendable to 5 or 6 levels with A 601 racks

EAN: 4002516050438 / Material number: 10973890 / Old Material Number: 69650101D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Mobile unit for inserts	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Capacity	
DIN mesh trays [number]	18
SPRI mesh trays [number]	12
Dimensions and weight	
External dimensions, net height in in. (mm)	26 (659)
External dimensions, net width in in. (mm)	25 5/8 (650)
External dimensions, net depth in in. (mm)	28 (710)
External dimensions, gross height in in. (mm)	31 5/16 (795)
External dimensions, gross width in in. (mm)	33 7/8 (860)
External dimensions, gross depth in in. (mm)	28 3/16 (715)
Net weight in lb (kg)	22 (9.80)
Gross weight in lb (kg)	69 (31.00)