

A 314

Mesh trays

For the optimum reprocessing of ophthalmic lumened instruments.



- Incl. 1 nozzle manifold to connect ophthalmic instruments
- Mesh tray with removable lid
- Particularly suitable for A 207 mobile unit

EAN: 4002516321262 / Material number: 10979670 / Old Material Number: 69631400D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Other inserts	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of ophthalmology instruments	•
Reprocessing of lumened instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	1 1/2 (38) in (mm)
External dimensions, net width in mm	8 13/16 (223) in (mm)
External dimensions, net depth in mm	6 (151) in (mm)
External dimensions, gross height in mm	3 (75)
External dimensions, gross width in mm	7 1/8 (180)
External dimensions, gross depth in in (mm)	9 11/16 (245)
Net weight in kg	2 (0.72)
Gross weight in kg	2 (0.85)