

A 505

Mobile unit

With 3 levels for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays

EAN: 4002516062998 / Material number: 10998970 / Old Material Number: 69650500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Mobile unit for inserts	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
Capacity	
DIN mesh trays [number]	9
SPRI mesh trays [number]	6
Dimensions and weight	
External dimensions, net height in in. (mm)	21 7/16 (543)
External dimensions, net width in in. (mm)	25 5/8 (650)
External dimensions, net depth in in. (mm)	31 1/8 (790)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	10 (4.40)
Gross weight in lb (kg)	15 (6.40)