


A 605 Module

For optimum reprocessing of vials and test tubes, etc.



- Equipped with 98 injector nozzles
- Max. H 160 mm, max. Ø 20 mm
- For use in A 503 mobile unit

EAN: 4002516063339 / Material number: 10999010 / Old Material Number: 69660500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert / module for test tubes	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of pipettes / phials	•
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	98
Vials [number]	98
Dimensions and weight	
External dimensions, net height in mm	6 1/4 (158) in (mm)
External dimensions, net width in mm	10 (254) in (mm)
External dimensions, net depth in mm	23 13/16 (604) in (mm)
External dimensions, gross height in mm	10 11/16 (270)
External dimensions, gross width in mm	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in kg	10 (4.40)
Gross weight in kg	15 (6.40)