

A 612

Module

For the optimum reprocessing of large-volume laboratory glassware.



- Can be equipped with injector nozzles based on individual needs
- Suitable for 2–10 l laboratory flasks, 2–6 l round flasks
- Suitable for 2–5 l Erlenmeyer flasks, 2–5 l measuring flasks
- For use in A 503 transfer cart

EAN: 4002516063353 / Material number: 10999050 / Old Material Number: 69661200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert / module for narrow necked glassware	•
Module for large-volume laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Large-volume laboratory glassware [number]	4
Dimensions and weight	
External dimensions, net height in mm	8 3/4 (221) in (mm)
External dimensions, net width in mm	10 (253) in (mm)
External dimensions, net depth in mm	23 15/16 (608) in (mm)
External dimensions, gross height in mm	9 11/16 (245)
External dimensions, gross width in mm	10 11/16 (270)
External dimensions, gross depth in in (mm)	25 (635)
Net weight in kg	8 (3.60)
Gross weight in kg	10 (4.40)