

A 506

Mobile unit

With 4 levels for cleaning laboratory glassware and utensils.



- For holding inserts and mesh trays

EAN: 4002516063223 / Material number: 10999130 / Old Material Number: 69650600D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Mobile unit for inserts	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
Capacity	
DIN mesh trays [number]	12
SPRI mesh trays [number]	8
Dimensions and weight	
External dimensions, net height in mm	21 7/8 (555) in (mm)
External dimensions, net width in mm	25 5/8 (650) in (mm)
External dimensions, net depth in mm	31 1/8 (790) in (mm)
External dimensions, gross height in mm	10 11/16 (270)
External dimensions, gross width in mm	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in kg	10 (4.40)
Gross weight in kg	15 (6.40)