


A 623

Module

For the optimum loading of measuring cylinders.



- Equipped with 4 nozzles
- For measuring cylinders with 1–2 l capacity
- Pivoting supports for ease of use

EAN: 4002516071440 / Material number: 11025560 / Old Material Number: 69662300D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert / module for measuring cylinder	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Capacity	
Measuring cylinders [number]	4
Dimensions and weight	
External dimensions, net height in in. (mm)	17 1/4 (437)
External dimensions, net width in in. (mm)	10 1/16 (255)
External dimensions, net depth in in. (mm)	23 15/16 (608)
External dimensions, gross height in in. (mm)	21 5/16 (540)
External dimensions, gross width in in. (mm)	20 7/8 (530)
External dimensions, gross depth in in. (mm)	19 11/16 (500)
Net weight in lb (kg)	14 (6.30)
Gross weight in lb (kg)	18 (7.80)