


A 847

Injector nozzle

With support, for cleaning Erlenmeyer flasks.



- For Erlenmeyer flasks, up to 5 L
- Nozzle length: 127 mm
- Nozzle Ø: 10 mm, Ø of support: 142 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516071471 / Material number: 11025600 / Old Material Number: 69684700D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	6 3/4 (170) in (mm)
External dimensions, net width in mm	9 13/16 (249) in (mm)
External dimensions, net depth in mm	1 7/16 (35) in (mm)
External dimensions, gross height in mm	12 1/16 (305)
External dimensions, gross width in mm	11 13/16 (300)
External dimensions, gross depth in in (mm)	7 7/8 (200)
Net weight in kg	1 (0.12)
Gross weight in kg	1 (0.27)