

A 848

Injector nozzle

With support, for cleaning measuring flasks.



- For measuring flasks, up to 5 L
- Nozzle length: 350 mm
- Nozzle Ø: 10 mm, Ø of support: 142 mm
- Suitable for A 313 and A 612 modules

EAN: 4002516071488 / Material number: 11025610 / Old Material Number: 69684800D

|  |                          |
|--|--------------------------|
| <b>Product category</b>                        |                          |
| Glassware Washers, laboratory                  | •                        |
| Large-capacity washer-disinfectors, laboratory | •                        |
| <b>Machine affiliation</b>                     |                          |
| PG 8504  | •                        |
| PG 8583  | •                        |
| PG 8583 CD                                     | •                        |
| PG 8593  | •                        |
| PLW 8615                                       | •                        |
| PLW 8616                                       | •                        |
| PLW 8617                                       | •                        |
| <b>Product category</b>                        |                          |
| Injector nozzle                                | •                        |
| <b>Product characteristics</b>                 |                          |
| Material                                       | Stainless steel, Plastic |
| Colour   | Stainless Steel, Grey    |
| <b>Application area</b>                        |                          |
| Reprocessing of laboratory glassware           | •                        |
| <b>Dimensions and weight</b>                   |                          |
| External dimensions, net height in mm          | 13 13/16 (350) in (mm)   |
| External dimensions, net width in mm           | 0 13/16 (20) in (mm)     |
| External dimensions, net depth in mm           | 0 13/16 (20) in (mm)     |
| External dimensions, gross height in mm        | 13 13/16 (350)           |
| External dimensions, gross width in mm         | 7 1/8 (180)              |
| External dimensions, gross depth in in (mm)    | 2 (50)                   |
| Net weight in kg                               | 1 (0.30)                 |
| Gross weight in kg                             | 2 (0.53)                 |