

A 306/1

Module

For the optimum loading of large-volume laboratory glassware.



- Suitable for 4 1–2 litre measuring cylinders, high shape
- Suitable for 2 1–2 litre measuring cyl., low shape + 2 high shape
- Plastic-coated contact surfaces
- Measuring cylinders even removable when module is in machine
- Maximum length of measuring cylinders: 500 mm
- Up to 2 modules can be inserted in baskets A 150 and APLW 106

EAN: 4002516076254 / Material number: 11035210 / Old Material Number: 69630601D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert / module for measuring cylinder	•
Module for large-volume laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Colour	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of measuring cylinders	•
Capacity	
Measuring cylinders [number]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	235
External dimensions, net width in mm	418
External dimensions, net depth in mm	471
External dimensions, gross height in mm	10 11/16 (270)
External dimensions, gross width in mm	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in kg	5 (2.10)
Gross weight in kg	7 (3.10)