

APH 592
 Colour transfer belt
 for the APH 510 PRT110 protocol printer.



EAN: 8809166675907 / Material number: 11053500 / Old Material Number: 68817005D

Product category	
Glassware Washers, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
PWD 8541	•
PWD 8545	•
PWD 8546	•
PWD 8549	•
Product category	
Process documentation / operating documentation / communication	•
Product characteristics	
Colour	Red, Black
Application area	
Reprocessing of pipettes / phials	•
Reprocessing of laboratory glassware	•
Reprocessing of slides	•
Reprocessing of butyrometers	•
Reprocessing of pharmacy / infusion bottles	•
Reprocessing anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of OR shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of orthopaedic surgery instruments	•
Reprocessing of intramedullary nails / bone marrow drills	•
Reprocessing of trays	•
Reprocessing of mesh trays	•
Reprocessing of measuring cylinders	•
Reprocessing of containers, e.g., soda lime containers	•
Reprocessing of bedpans and urine bottles	•

APH 592
 Colour transfer belt
 for the APH 510 PRT110 protocol printer.

EAN: 8809166675907 / Material number: 11053500 / Old Material Number: 68817005D



Reprocessing of different utensils

•

Dimensions and weight

External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	2 1/4 (55)
External dimensions, gross width in mm	7 1/8 (180)
External dimensions, gross depth in in (mm)	9 7/8 (250)
Net weight in kg	1 (0.04)
Gross weight in kg	1 (0.05)

Standards / certs / symbols

CE

•