

A 312
Module
For reprocessing viscometers.



- Suitable for up to 8.0 viscometers
- Equipped with 8 x A 867 special nozzles
- Nozzle spacing 118 mm across the width, 118 mm across the depth

EAN: 4002516087397 / Material number: 11054470 / Old Material Number: 69631200D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Module for viscometers	•
Product characteristics	
Material	Stainless steel, Plastic, Silicone
Colour	Stainless Steel, Black
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Viscometer [quantity]	8.0
Dimensions and weight	
External dimensions, net height in mm	2 5/8 (66) in (mm)
External dimensions, net width in mm	5 1/4 (133) in (mm)
External dimensions, net depth in mm	18 7/8 (479) in (mm)
External dimensions, gross height in mm	9 7/8 (250)
External dimensions, gross width in mm	9 1/16 (230)
External dimensions, gross depth in in (mm)	19 15/16 (505)
Net weight in kg	4 (1.58)
Gross weight in kg	6 (2.33)