


A 836

Set of injector nozzles

For optimum reprocessing of hollow instruments used in dental applications.



- For use in A 105/1 upper basket
- Packaging unit with 22 dental injector nozzles (A 833)

EAN: 4002516204091 / Material number: 11054660 / Old Material Number: 69783600D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	4 3/4 (120)
External dimensions, gross depth in in (mm)	6 3/4 (170)
Net weight in kg	1 (0.20)
Gross weight in kg	1 (0.20)