


A 837

Irrigation sleeve set

For reprocessing hollow instruments used in ENT applications in particular.



- For use in A 105/1 upper basket
- Incl. cap and mounting clip
- Packaging unit with 22 irrigation sleeves (A 817)

EAN: 4002516204107 / Material number: 11054690 / Old Material Number: 69783700D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Irrigation sleeve	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Application area</b>	
Reprocessing of ENT instruments	•
Reprocessing of lumened instruments	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	2 3/16 (55)
External dimensions, gross width in mm	7 1/8 (180)
External dimensions, gross depth in in (mm)	9 7/8 (250)
Net weight in kg	2 (0.53)
Gross weight in kg	2 (0.53)