

A 315 Module

For increased capacity and reprocessing of hollow instruments.



- Increases the injector capacity from 22 to 44
- Backwards-compatible with all A 151 products on the market
- Compatible with all standard dental & ENT adapters & accessories

EAN: 4002516087861 / Material number: 11055980 / Old Material Number: 69631500D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert for ENT instruments	•
Insert / module for dental instruments	•
Other inserts / other modules	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Option to use a central filter	•
Application area	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of lumened instruments	•
Capacity	
ENT instruments [number]	22
Handpieces [number]	22
Siphons [number]	22
Lumened instruments [number]	22
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	145
External dimensions, net width in mm	492
External dimensions, net depth in mm	542
External dimensions, gross height in mm	10 11/16 (270)
External dimensions, gross width in mm	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in kg	3 (1.35)
Gross weight in kg	6 (2.35)
Accessory pack	
A 838 key set [number]	1