

A 841

Injector nozzle

With plastic supports, nozzle diameter 6, length 210 mm, 1 item



- Suitable for injector carts and baskets, with screw connection
- Part of SimpleLoad system
- Packaging unit with 1 nozzle
- Ideal for load items with neck diameters of between 12 mm and 85 mm
- Diameter of cap at top: 11 mm

EAN: 4002516188759 / Material number: 11056190 / Old Material Number: 69684100D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	8 5/16 (210) in (mm)
External dimensions, net width in mm	3 5/8 (92) in (mm)
External dimensions, net depth in mm	3 7/16 (86) in (mm)
External dimensions, gross height in mm	4 1/4 (107)
External dimensions, gross width in mm	5 1/8 (130)
External dimensions, gross depth in in (mm)	7 11/16 (195)
Net weight in kg	1 (0.03)
Gross weight in kg	1 (0.14)