

A 844

Injector nozzle

With plastic supports, nozzle diameter 2.5, length 80 mm, 1 item.



- Suitable for injector carts and baskets, with screw connection
- Part of SimpleLoad system
- Packaging unit with 1 nozzle

EAN: 4002516188780 / Material number: 11056600 / Old Material Number: 69684400D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Injector nozzle	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in in. (mm)	3 3/16 (80)
External dimensions, net width in in. (mm)	2 3/8 (60)
External dimensions, net depth in in. (mm)	2 3/16 (55)
External dimensions, gross height in in. (mm)	1 5/8 (40)
External dimensions, gross width in in. (mm)	2 9/16 (65)
External dimensions, gross depth in in (mm)	4 11/16 (118)
Net weight in lb (kg)	1 (0.01)
Gross weight in lb (kg)	1 (0.03)