

A 300/2

Module

For optimum loading of laboratory glassware – can be equipped with 8 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 8 injector nozzles
- Part of SimpleLoad system
- Nozzle spacing 118 mm across the width, 118 mm across the depth
- Particularly suitable for A 840 and A 841 (nozzle diameter 6 mm)

EAN: 4002516087632 / Material number: 11056770 / Old Material Number: 69630002D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert / module for narrow necked glassware	•
Insert for various bottles	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel, Black
Part of SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	8
Standard electrical connection	
Number of phases	0
Dimensions and weight	
External dimensions, net height in in. (mm)	2 7/8 (72)
External dimensions, net width in in. (mm)	5 1/4 (133)
External dimensions, net depth in in. (mm)	18 7/8 (479)
External dimensions, gross height in in. (mm)	7 7/8 (200)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	2 (0.56)
Gross weight in lb (kg)	4 (1.64)