


A 302/2 Module

For optimum loading of laboratory glassware – can be equipped with 32 nozzles.



- Suitable for lab. glassware such as Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 32 injector nozzles
- Part of SimpleLoad system
- Nozzle spacing 60 mm across the width, 57 mm across the depth
- Particularly suitable for A 844 and A 845 (nozzle diameter 2.5 mm)

EAN: 4002516087656 / Material number: 11056790 / Old Material Number: 69630202D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert / module for narrow necked glassware	•
Insert for various bottles	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	32
Standard electrical connection	
Number of phases	0
Dimensions and weight	
External dimensions, net height in in. (mm)	2 7/8 (72)
External dimensions, net width in in. (mm)	7 11/16 (195)
External dimensions, net depth in in. (mm)	18 7/8 (479)
External dimensions, gross height in in. (mm)	7 7/8 (200)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	3 (1.09)
Gross weight in lb (kg)	5 (2.16)