


A 862

Holding frames

For injector modules A 302/2 and A 302/3



- Arrangement: 4 x 8 fields
- For use on A 302/2 and A 302/3 injector modules
- For the optimum positioning of laboratory glassware
- Plastic-coated

EAN: 4002516087687 / Material number: 11056820 / Old Material Number: 69686200D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Other inserts / other modules	•
Product characteristics	
Material	Stainless steel, Plastic, Wire, rilsan-coated
Colour	Stainless Steel, Black
Number of compartments	32
Compartment width in mm	2 1/4 (56)
Compartment length in mm	2 1/16 (52)
Part of SimpleLoad system	•
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	1 3/16 (30)
External dimensions, net width in in. (mm)	9 7/8 (250)
External dimensions, net depth in in. (mm)	13 13/16 (350)
External dimensions, gross height in in. (mm)	1 7/16 (35)
External dimensions, gross width in in. (mm)	11 13/16 (300)
External dimensions, gross depth in in. (mm)	15 3/4 (400)
Net weight in lb (kg)	1 (0.10)
Gross weight in lb (kg)	1 (0.18)