

APLW 052/1

Set of injector nozzles

With plastic supports, Ø 4 mm, length 175 mm.

EAN: 4002516130918 / Material number: 11094000 / Old Material Number: 69305201D



- Suitable for e.g., lower basket APLW 095
- Packaging unit with 39 injector nozzles
- Suitable e.g., for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 5 mm

<b>Product category</b>	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
PLW 8505	•
<b>Product category</b>	
Injector nozzle	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	7 1/2 (190)
External dimensions, gross width in mm	12 1/4 (310)
External dimensions, gross depth in in (mm)	12 1/4 (310)
Net weight in kg	2 (0.78)
Gross weight in kg	3 (1.05)