


APLW 002

Mobile unit

For optimum loading of large-volume laboratory glassware, 5 and 10 l.



- Equipped with 5 nozzles
- For loading up to 4 bottles with a max. Ø of 240 mm
- For loading up to 5 bottles with a max. Ø of 190 mm

EAN: 4002516115304 / Material number: 11104490 / Old Material Number: 69300200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Module for large-volume laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Flasks [number]	5
Dimensions and weight	
External dimensions, net height in mm	413
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	22 7/8 (580)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	14 (5.94)
Gross weight in kg	19 (8.49)