

APLW 004

Mobile unit

For optimum loading of large-volume laboratory glassware, up to 50 l.



- Equipped with 1 nozzle

EAN: 4002516115427 / Material number: 11104510 / Old Material Number: 69300400D

| | |
|--|--------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6111 | • |
| Product category | |
| Injector basket | • |
| Module for large-volume laboratory glassware | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Flasks [number] | 1 |
| Dimensions and weight | |
| External dimensions, net height in mm | 588 |
| External dimensions, net width in mm | 535 |
| External dimensions, net depth in mm | 582 |
| External dimensions, gross height in mm | 30 (760) |
| External dimensions, gross width in mm | 24 5/8 (625) |
| External dimensions, gross depth in in (mm) | 28 1/4 (715) |
| Net weight in kg | 11 (4.68) |
| Gross weight in kg | 22 (9.83) |