

APLW 006

Mobile unit

For optimum loading of up to 121 pipettes.



- For loading a max. of 121 pipettes in a vertical arrangement

EAN: 4002516115441 / Material number: 11104530 / Old Material Number: 69300600D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Application area	
Reprocessing of pipettes / phials	•
Capacity	
Pipettes [number]	121
Dimensions and weight	
External dimensions, net height in mm	195
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/8 (610)
Net weight in kg	11 (4.94)
Gross weight in kg	17 (7.38)