


APLW 008

Mobile unit

For optimum loading of up to 3 pipette cartridges.



- For lower level
- Size of pipette cartridge: W 84 x H 384 x D 84 mm

EAN: 4002516115465 / Material number: 11104550 / Old Material Number: 69300800D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6111	•
Product category	
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Application area	
Reprocessing of pipettes / phials	•
Dimensions and weight	
External dimensions, net height in mm	456
External dimensions, net width in mm	179
External dimensions, net depth in mm	366
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/8 (610)
Net weight in kg	17 (7.64)
Gross weight in kg	22 (9.64)