

APLW 010

Upper basket

With 156 unassembled injector positions for nozzles, Ø 2.5 mm.



**Image  
in progress**

- Only suited for use on the top levels

EAN: 4002516115489 / Material number: 11104570 / Old Material Number: 69301000D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow necked glassware [number]	156
<b>Dimensions and weight</b>	
External dimensions, net height in mm	162
External dimensions, net width in mm	534
External dimensions, net depth in mm	569
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	12 (5.16)
Gross weight in kg	17 (7.58)