

APLW 012

Upper basket

With 84 unassembled injector positions for nozzles, Ø 2.5 mm.



- Only suited for use on the top levels

EAN: 4002516115502 / Material number: 11104590 / Old Material Number: 69301200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	84
Dimensions and weight	
External dimensions, net height in mm	158
External dimensions, net width in mm	532
External dimensions, net depth in mm	570
External dimensions, gross height in mm	13 7/8 (350)
External dimensions, gross width in mm	22 1/8 (560)
External dimensions, gross depth in in (mm)	23 5/8 (600)
Net weight in kg	11 (4.60)
Gross weight in kg	17 (7.40)