

APLW 013  
Upper basket  
With 42 unassembled injector positions.



- Only suited for use on the top levels

EAN: 4002516115519 / Material number: 11104600 / Old Material Number: 69301300D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow necked glassware [number]	42
<b>Dimensions and weight</b>	
External dimensions, net height in mm	160
External dimensions, net width in mm	534
External dimensions, net depth in mm	569
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	11 (4.55)
Gross weight in kg	16 (6.95)