

APLW 013/2

Upper basket

Assembled with 32 injector nozzles APLW 051 and 10 injector nozzles APLW 055



**Image  
in progress**

- Only suited for use on the top levels

EAN: 4002516115533 / Material number: 11104620 / Old Material Number: 69301302D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow necked glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
<b>Dimensions and weight</b>	
External dimensions, net height in mm	250
External dimensions, net width in mm	535
External dimensions, net depth in mm	570
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	11 (4.55)
Gross weight in kg	20 (9.00)
<b>Accessory pack</b>	
APLW 051, injector nozzle Ø 4 mm, length 110 mm [number]	32
APLW 055, injector nozzle Ø 6 mm, length 175 mm [number]	10