

APLW 014  
 Upper basket  
 With 20 unassembled injector positions.



- Only suited for use on the top levels

EAN: 4002516115540 / Material number: 11104630 / Old Material Number: 69301400D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow necked glassware [number]	20
Flasks [number]	20
Erlenmeyer flasks [number]	20
Round flasks [number]	20
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0 7/8 (21) in (mm)
External dimensions, net width in mm	0 13/16 (20) in (mm)
External dimensions, net depth in mm	0 11/16 (16) in (mm)
External dimensions, gross height in mm	13 7/16 (340)
External dimensions, gross width in mm	21 7/8 (555)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	10 (4.22)
Gross weight in kg	15 (6.65)