

APLW 015
Upper basket
With 16 unassembled injector positions.



- Only suited for use on the top levels

EAN: 4002516115571 / Material number: 11104650 / Old Material Number: 69301500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Capacity	
Narrow necked glassware [number]	16
Flasks [number]	16
Erlenmeyer flasks [number]	16
Round flasks [number]	16
Dimensions and weight	
External dimensions, net height in in. (mm)	6 5/16 (160)
External dimensions, net width in in. (mm)	21 (532)
External dimensions, net depth in in. (mm)	22 1/16 (560)
External dimensions, gross height in in. (mm)	13 7/16 (340)
External dimensions, gross width in in. (mm)	22 1/16 (560)
External dimensions, gross depth in in. (mm)	23 5/8 (600)
Net weight in lb (kg)	9 (4.00)
Gross weight in lb (kg)	15 (6.60)