

APLW 015
Upper basket
With 16 unassembled injector positions.



- Only suited for use on the top levels

EAN: 4002516115571 / Material number: 11104650 / Old Material Number: 69301500D

| | |
|--|------------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector basket | • |
| Upper basket | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Colour | Blue, Stainless Steel, White |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Narrow necked glassware [number] | 16 |
| Flasks [number] | 16 |
| Erlenmeyer flasks [number] | 16 |
| Round flasks [number] | 16 |
| Dimensions and weight | |
| External dimensions, net height in mm | 160 |
| External dimensions, net width in mm | 532 |
| External dimensions, net depth in mm | 560 |
| External dimensions, gross height in mm | 13 1/2 (340) |
| External dimensions, gross width in mm | 22 1/8 (560) |
| External dimensions, gross depth in in (mm) | 23 5/8 (600) |
| Net weight in kg | 9 (4.00) |
| Gross weight in kg | 15 (6.60) |