

APLW 017
Upper basket
With 121 unassembled injector positions.



- Only suited for use on the top levels
- 200 x 490 mm chamber

EAN: 4002516115595 / Material number: 11104670 / Old Material Number: 69301700D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Vials [number]	121
Dimensions and weight	
External dimensions, net height in mm	160
External dimensions, net width in mm	532
External dimensions, net depth in mm	560
External dimensions, gross height in mm	13 7/8 (350)
External dimensions, gross width in mm	22 1/8 (560)
External dimensions, gross depth in in (mm)	24 1/8 (610)
Net weight in kg	10 (4.20)
Gross weight in kg	16 (7.00)