

APLW 018
 Upper basket
 With 121 + 24 unassembled injector positions.



- Only suited for use on the top levels

EAN: 4002516115601 / Material number: 11104680 / Old Material Number: 69301800D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Vials [number]	121
Narrow necked glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
Dimensions and weight	
External dimensions, net height in mm	162
External dimensions, net width in mm	534
External dimensions, net depth in mm	569
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	12 (5.33)
Gross weight in kg	18 (7.76)