


APLW 019/1

Upper basket

Assembled with 121 injector nozzles APLW 047.



- Only suited for use on the top levels

EAN: 4002516115731 / Material number: 11104710 / Old Material Number: 69301901D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Vials [number]	121
Dimensions and weight	
External dimensions, net height in mm	200
External dimensions, net width in mm	535
External dimensions, net depth in mm	570
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	15 (6.64)
Gross weight in kg	20 (9.06)
Accessory pack	
APLW 047 injector nozzle, Ø 2.5 mm, length 50 mm [number]	121