

APLW 020/1
 Lower basket
 equipped with 121 APLW 047 injector nozzles.



- Only suited for use on the lowest level

EAN: 4002516115755 / Material number: 11104730 / Old Material Number: 69302001D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	121
Vials [number]	121
Dimensions and weight	
External dimensions, net height in mm	200
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (562)
External dimensions, gross depth in in (mm)	26 5/8 (675)
Net weight in kg	20 (8.64)
Gross weight in kg	25 (11.26)
Accessory pack	
APLW 047 injector nozzle, Ø 2.5 mm, length 50 mm [number]	121