

APLW 021

Lower basket

With 110 unassembled injector positions for nozzles, Ø 2.5 mm.



**Image
in progress**

- Only suited for use on the lowest level

EAN: 4002516115762 / Material number: 11104740 / Old Material Number: 69302100D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	110
Vials [number]	110
Narrow necked glassware [number]	110
Flasks [number]	110
Dimensions and weight	
External dimensions, net height in mm	575
External dimensions, net width in mm	504
External dimensions, net depth in mm	450
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (562)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	14 (5.91)
Gross weight in kg	19 (8.35)