

APLW 023

Lower basket

With 56 unassembled injector positions for nozzles, Ø 2.5 mm.



**Image
in progress**

- Only suited for use on the lowest level

EAN: 4002516115786 / Material number: 11104760 / Old Material Number: 69302300D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	56
Flasks [number]	56
Erlenmeyer flasks [number]	56
Round flasks [number]	56
Dimensions and weight	
External dimensions, net height in mm	575
External dimensions, net width in mm	504
External dimensions, net depth in mm	480
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (562)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	13 (5.48)
Gross weight in kg	18 (7.91)