

APLW 024/2

Lower basket

equipped with 32 APLW 051 injector nozzles and 10 APLW 055 injector nozzles.



**Image
in progress**

- Only suited for use on the lowest level

EAN: 4002516115816 / Material number: 11104790 / Old Material Number: 69302402D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
Dimensions and weight	
External dimensions, net height in mm	220
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	13 7/8 (352)
External dimensions, gross width in mm	22 1/4 (562)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	14 (5.97)
Gross weight in kg	24 (10.48)
Accessory pack	
APLW 051, injector nozzle Ø 4 mm, length 110 mm [number]	32
APLW 055, injector nozzle Ø 6 mm, length 175 mm [number]	10