

APLW 024/3

Lower basket

Assembled with 30 injector nozzles APLW 066 and 12 injector nozzles APLW 067.



**Image  
in progress**

- Only suited for use on the lowest level

EAN: 4002516115823 / Material number: 11104800 / Old Material Number: 69302403D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow necked glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
<b>Dimensions and weight</b>	
External dimensions, net height in mm	220
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	13 7/8 (352)
External dimensions, gross width in mm	22 1/4 (562)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	12 (5.22)
Gross weight in kg	31 (13.87)
<b>Accessory pack</b>	
APLW 066 injector nozzle, Ø 6 mm, length 175 mm [number]	30
APLW 067 injector nozzle, Ø 6 mm, length 225 mm [number]	12