

APLW 024/3

Lower basket

Assembled with 30 injector nozzles APLW 066 and 12 injector nozzles APLW 067.



- Only suited for use on the lowest level

EAN: 4002516115823 / Material number: 11104800 / Old Material Number: 69302403D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
Capacity	
Narrow necked glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
Dimensions and weight	
External dimensions, net height in in. (mm)	8 11/16 (220)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	22 11/16 (575)
External dimensions, gross height in in. (mm)	13 7/8 (352)
External dimensions, gross width in in. (mm)	22 3/16 (562)
External dimensions, gross depth in in. (mm)	24 1/4 (615)
Net weight in lb (kg)	12 (5.22)
Gross weight in lb (kg)	31 (13.88)
Accessory pack	
APLW 066 injector nozzle, Ø 6 mm, length 175 mm [number]	30
APLW 067 injector nozzle, Ø 6 mm, length 225 mm [number]	12