

APLW 025

Lower basket

With 24 unassembled injector positions for nozzles.



**Image**  
in progress

- Only suited for use on the lowest level
- 230 x 490 mm chamber

EAN: 4002516115830 / Material number: 11104810 / Old Material Number: 69302500D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	24
Flasks [number]	24
Erlenmeyer flasks [number]	24
Round flasks [number]	24
Dimensions and weight	
External dimensions, net height in mm	136
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	14 1/4 (360)
External dimensions, gross width in mm	22 1/8 (560)
External dimensions, gross depth in in (mm)	23 5/8 (600)
Net weight in kg	10 (4.40)
Gross weight in kg	17 (7.40)