


APLW 026/1

Lower basket

Assembled with 20 injector nozzles APLW 056.

EAN: 4002516115861 / Material number: 11104840 / Old Material Number: 69302601D



- Only suited for use on the lowest level

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	20
Flasks [number]	20
Erlenmeyer flasks [number]	20
Round flasks [number]	20
Dimensions and weight	
External dimensions, net height in mm	270
External dimensions, net width in mm	535
External dimensions, net depth in mm	570
External dimensions, gross height in mm	13 1/2 (340)
External dimensions, gross width in mm	21 7/8 (555)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	11 (4.92)
Gross weight in kg	21 (9.20)
Accessory pack	
APLW 056 injector nozzle, Ø 6 mm, length 225 mm [number]	20