

APLW 027

Lower basket

With 16 unassembled injector positions for nozzles.



**Image
in progress**

- Only suited for use on the lowest level

EAN: 4002516115878 / Material number: 11104850 / Old Material Number: 69302700D

| | |
|--|-------------------------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector basket | • |
| Lower basket | • |
| Product characteristics | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black, White |
| Application area | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Narrow necked glassware [number] | 16 |
| Flasks [number] | 16 |
| Erlenmeyer flasks [number] | 16 |
| Round flasks [number] | 16 |
| Dimensions and weight | |
| External dimensions, net height in mm | 570 |
| External dimensions, net width in mm | 536 |
| External dimensions, net depth in mm | 516 |
| External dimensions, gross height in mm | 13 1/2 (340) |
| External dimensions, gross width in mm | 21 7/8 (555) |
| External dimensions, gross depth in in (mm) | 24 1/4 (615) |
| Net weight in kg | 11 (4.74) |
| Gross weight in kg | 16 (7.05) |