

APLW 029

Lower basket

With 27 unassembled injector positions for nozzles.



- Only suited for use on the lowest level

EAN: 4002516115892 / Material number: 11104870 / Old Material Number: 69302900D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow necked glassware [number]	12
Flasks [number]	12
Erlenmeyer flasks [number]	12
Round flasks [number]	12
<b>Dimensions and weight</b>	
External dimensions, net height in mm	137
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	11 (4.74)
Gross weight in kg	16 (7.18)