

APLW 030

Lower basket

With 12 unassembled injector positions for nozzles.



**Image
in progress**

- Only suited for use on the lowest level
- 220 x 490 mm chamber

EAN: 4002516115908 / Material number: 11104880 / Old Material Number: 69303000D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Application area	
Reprocessing of laboratory glassware	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black, White
Capacity	
Narrow necked glassware [number]	12
Flasks [number]	12
Erlenmeyer flasks [number]	12
Round flasks [number]	12
Dimensions and weight	
External dimensions, net height in in. (mm)	5 3/8 (135)
External dimensions, net width in in. (mm)	16 3/4 (424)
External dimensions, net depth in in. (mm)	22 15/16 (582)
External dimensions, gross height in in. (mm)	14 (355)
External dimensions, gross width in in. (mm)	22 1/4 (565)
External dimensions, gross depth in in. (mm)	24 1/4 (615)
Net weight in lb (kg)	10 (4.32)
Gross weight in lb (kg)	15 (6.72)