


APLW 030/1

Lower basket

equipped with 6 APLW 067 injector nozzles and 6 APLW 066 injector nozzles.

EAN: 4002516115915 / Material number: 11104890 / Old Material Number: 69303001D



- Only suited for use on the lowest level

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	12
Flasks [number]	12
Erlenmeyer flasks [number]	12
Round flasks [number]	12
Dimensions and weight	
External dimensions, net height in mm	235
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	14 (355)
External dimensions, gross width in mm	22 1/4 (565)
External dimensions, gross depth in in (mm)	24 1/4 (615)
Net weight in kg	12 (5.33)
Gross weight in kg	19 (8.55)
Accessory pack	
APLW 066 injector nozzle, Ø 6 mm, length 175 mm [number]	6
APLW 067 injector nozzle, Ø 6 mm, length 225 mm [number]	6